

Mouse IL17RB Recombinant Protein



RPPB0619

Product Information Protein Information

Product SKU:

RPPB0619

Accession:

Q9Z132

Host:

Sf9, Baculovirus cells.

Protein description:

IL17RB produced in Sf9 Baculovirus cells is a single, glycosylated polypeptide chain containing 511 amino acids (18-286 a.a.) and having a molecular mass of 57.1kDa (Molecular size on SDS-PAGE will appear at approximately 57-70kDa). IL17RB is expressed with an 242 amino acid hlgG-His tag at C-Terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered colorless solution.

Synonyms:

IL17rb, Evi27, IL-17ER, IL-17Rh1, IL17br, IL17RH1, IL-17 receptor B, IL-17RB, Interleukin-17 receptor B.

Formulation:

IL17RB protein solution (1mg/ml) contains Phosphate Buffered Saline (pH 7.4) and 10% glycerol.

Purity:

Greater than 90.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Amino Acid Sequence:

ADPREPTIQC GSETGPSPEW MVQHTLTPGD LRDLQVELVK TSVAEEFSI LMNISWILRA DASIRLLKAT
KICVSGKNNM NSYSCVRCNY TEAFQSQTRP SGGKWTFYSYV GFPVELSTLY LISAHNIPNA NMNEDSPSL
VNFTSPGCLN HVMKYKKQCT EAGSLWDPDI TACKKNEKMV EVNFTTNPLG NRYTILIQRD TTLGFSRVLE
NKL MRTSVAI PVTEESEGA VQLTPYLHTC GND CIRREGT VVLCSETSAP IPPDDNRRML GGVEPKSCDK
THTCPPCPAP ELLGGPSVFL FPPKPKDTLM ISRTPEVTCV VVDVSHEDPE VKFNWYVDGV EVHNAKTKPR
EEQYNSTYRV VSVLTVLHQD WLNGKEYKCK VSNKALPAPI EKTISKAKGQ PREPQVYTLPSRDELTKNQ
VSLTCLVKGF YPSDIAVEWE SNGQPENNYK TTPVLDSDG SFFLYSKLTV DKSRWQQGNV FSCSV MHEAL
HNHYTQKSLS LSPGKHHHHH H