

# Human CD4 (203-317) Recombinant Protein



RPPB0890

---

## Product Information    Protein Information

**Product SKU:**

RPPB0890

**Protein description:**

CD-4 Human Recombinant His Tag protein encoding 203-317 amino acids.

**Host:**

Escherichia Coli.

**Appearance:**

Sterile lyophilized powder.

**Synonyms:**

gp55, HLA-2, L3 / T4, Ly-4, T cell antigen T4/LEU3, T4, sCD4.

**Formulation:**

Each mg contains 1X PBS, 0.1% SDS and 1mM DTT.

**Solubility:**

It is recommended to reconstitute the lyophilized CD-4 (203-317) in sterile 18MΩ-cm H<sub>2</sub>O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

**Stability:**

Store vial at -20oC to -80oC. When stored at the recommended temperature, this protein is stable for 12 months. Please avoid freeze-thaw cycles.