

Human TGFB3 (207 a.a.) Recombinant Protein



RPPB0979

Product Information Protein Information

Product SKU:

RPPB0979

Accession:

P10600

Host:

Escherichia Coli.

Protein description:

TGFB3 Human Recombinant protein encoding amino acids 644-850 and total molecular mass of 50 kDa (Including GST Tag).

Appearance:

Sterile Filtered clear solution.

Synonyms:

Transforming Growth Factor-beta3, TGFB3, ARVD, FLJ16571, TGF-beta3.

Formulation:

100µl purified human TGF Beta 3protein at 100µg/ml in 50mM Tris-Acetate, pH-7.5, 1mM EDTA and 20% Glycerol.

Stability:

TGF-beta 3 Human Recombinant although stable at 4°C for 1 week, should be stored at -20°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Please avoid freeze-thaw cycles.