

# Human Inhibin alpha A Chain Recombinant Protein



RPPB1284

---

## Product Information    Protein Information

**Product SKU:**

RPPB1284

**Host:**

Escherichia Coli.

**Protein description:**

Inhibin-Alpha A chain Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 134 amino acids fragment (233-366) and having an amino-terminal hexahistidine tag, having a total molecular weight of 19.2 kDa. The Inhibin-Alpha A chain is purified by standard chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Inhibin-Alpha, A-Inhibin Subunit, Inhibin Alpha Chain, INHA.

**Formulation:**

Inhibin-A alpha chain is supplied in 20mM Tris HCl (pH-8), 5mM EDTA and 50% glycerol.

**Purity:**

Greater than 95.0% as determined by:(a) Analysis by RP-HPLC.(b) Analysis by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.