

Human ALPP Recombinant Protein



RPPB1430

Product Information Protein Information

Product SKU:

RPPB1430

Accession:

P05187

Host:

Escherichia Coli.

Protein description:

Placental Alkaline Phosphatase Human Recombinant encoding 154-287 amino acids expressed in E.coli, shows a 41kDa band on SDS-PAGE (including GST tag). PLAP Human Recombinant is purified by proprietary chromatographic techniques.

Appearance:

Sterile filtered liquid.

Synonyms:

ALP, PLAP, Alkaline phosphatase placental type, EC 3.1.3.1, PLAP-1, Alkaline phosphatase Regan isozyme.

Formulation:

PLAP in 50mM Tris-Acetate, pH7.5, 1mM EDTA and 20% Glycerol.

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.