

Pyrococcus Fruriosus DUT Recombinant Protein



RPPB1603

Product Information Protein Information

Product SKU:

RPPB1603

Appearance:

Sterile filtered liquid formulation.

Host:

Escherichia Coli.

Synonyms:

Thermostable dUTPase, dUTPase.

Formulation:

dUTPase is supplied in 20mM Tris-HCl (pH 8.2), 1mM DTT, 0.1mM EDTA, 100mM KCl, 0.1% Nonidet P40, 0.1% Tween 20 and 50% glycerol at a concentration of 1000U/μl of the enzyme.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.