

Yeast GLRX1 Recombinant Protein



RPPB1699

Product Information Protein Information

Product SKU:

RPPB1699

Accession:

P25373

Host:

Escherichia Coli.

Protein description:

Glutaredoxin Saccharomyces cerevisiae Recombinant containing 6x His tag at C-Terminus produced in E.Coli is a single, non-glycosylated, Polypeptide chain having a molecular mass of 16 kDa.

Appearance:

Sterile Filtered clear colorless solution.

Synonyms:

Thioltransferase, GRX, GLRX1, GRX1, GRX-1, GLRX-1, Glutathione-dependent oxidoreductase 1, Glutaredoxin 1.

Formulation:

Glutaredoxin solution contains 25mM Tris-HCl pH-7.5 & 0.01% Na Azide.

Purity:

Greater than 90% as determined by SDS-PAGE.

Stability:

1 week at 2-10°C. For long term store at -20 to -80°C.