

# ENO2 Protein Recombinant Protein



RPPB1702

---

## Product Information    Protein Information

**Product SKU:**

RPPB1702

**Protein description:**

Human Neurone Specific Enolase produced in Human CNS having a molecular mass of 45kDa.

**Accession:**

P09104

**Appearance:**

Sterile Filtered clear solution.

**Host:**

Human CNS.

**Synonyms:**

Gamma-enolase, EC 4.2.1.11, 2-phospho-D-glycerate hydro-lyase, Neural enolase, Neuron-specific enolase, NSE, Enolase 2, ENO2.

**Formulation:**

The protein solution is in 0.01M NaH<sub>2</sub>PO<sub>4</sub> buffer pH 7.4 containing 0.15M NaCl and 0.005M MgSO<sub>4</sub>.

**Purity:**

Greater than 96.0%.

**Stability:**

Human NSE although stable at 4°C for 1 week, should be stored at -15°C.