

Bovine Enterokinase His Recombinant Protein



RPPB1710

Product Information Protein Information

Product SKU:

RPPB1710

Accession:

P98072

Host:

Escherichia Coli.

Protein description:

Enterokinase Bovine Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 241 amino acids with a 6 × His at C-terminus and having a molecular mass of 28.0kDa. The Enterokinase Bovine is purified by proprietary chromatographic techniques.

Appearance:

Sterile liquid solution.

Synonyms:

Enteropeptidase, EC 3.4.21.9, Enterokinase, Serine protease 7, ENTK, MGC133046.

Formulation:

Bovine EK is supplied in 50mM Tris-HCl, pH 8.0, 0.5M NaCl and 50% glycerol.

Stability:

One year when stored at -20°C. Please avoid freeze-thaw cycles.