

HPSE WB Recombinant Protein



RPPB1820

Product Information Protein Information

Product SKU:

RPPB1820

Accession:

Q9Y251

Protein description:

Recombinant Heparanase protein HPA1 is produced in CHO cells. The protein is purified by several orthogonal chromatography steps.

Formulation:

Concentration: 1 µg /ml
Content: 100 ng
Buffer: LDS-PAGE buffer [140 mM Tris buffer pH 8.5, 10% Glycerol, 2% LDS, 0.015% EDTA, 1.88% (v/v) of 1% Serva Blue G250 and 0.625% (v/v) of 1% Phenol red].