

Human HTATIP2 Recombinant Protein



RPPB1836

Product Information Protein Information

Product SKU:

RPPB1836

Accession:

Q9BUP3

Host:

Escherichia Coli.

Protein description:

HTATIP2 Recombinant Human produced in E.Coli is a single, non-glycosylated polypeptide chain containing 262 amino acids (1-242 a.a.) and having a molecular mass of 29.3 kDa. The HTATIP2 is fused to a 20 amino acid His-Tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile filtered colorless solution.

Synonyms:

TIP30, CC3, SDR44U1, HTATIP2, EC=1.1.1.-, HIV-1 TAT-interactive protein 2, FLJ26963.

Formulation:

HTATIP2 protein solution (1mg/ml) containing 20mM Tris-HCl pH-8, 1mM DTT & 10% glycerol.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Amino Acid Sequence:

MGSSHHHHHH SGLVPRGSH MAETEALSKL REDFRMQNKS VFILGASGET GRVLLKEILE QGLFSKVTLI
GRRKLTDFDEE AYKNVNQEVV DFEKLDDYAS AFQGHVDGFC CLGTTRGKAG AEGFVRVDRD YVLKSAELAK
AGGCKHFNLL SSKGADKSSN FLYLQVKGEV EAKVEELKFD RYSVFRPGVL LCDRQESRPG EWLVRKFFGS
LPDSWARGHS VPVTVVVRAM LNNVVRPRDK QMELLENKAI HDLGKAHGSL KP.