

Human PGAM1 Recombinant Protein



RPPB2093

Product Information Protein Information

Product SKU:

RPPB2093

Accession:

P18669

Host:

Escherichia Coli.

Protein description:

PGAM1 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 274 amino acids (1-254 a.a.) and having a molecular mass of 30.9 kDa. The PGAM1 is fused to a 20 amino acid His Tag at N-Terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear colorless solution.

Synonyms:

Phosphoglycerate mutase isozyme B, PGAM-B, PGAMA.

Formulation:

The PGAM1 1mg/ml protein solution contains 20mM Tris-HCl pH-8, 1mM DTT, and 10% glycerol.

Purity:

Greater than 90% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Amino Acid Sequence:

MGSSHHHHHH SGLVPRGSH MAAYKLVLR HGESAWNLEN RFGWYDADL SPAGHEEAKR GGQALRDAGY
EFDICFTSVQ KRAIRTLWTV LDAIDQMWLP VVRTWRLNER HYGGLTGLNK AETAAKHGEA QVKIWRRSYD
VPPPPMEPDH PFYSNISKDR RYADLTEDQL PSCESLKDTI ARALPFWNEE IVPQIKEGKR VLIAAHGNSL
RGIVKHLEGL SEEAIMELNL PTGIPIVYEL DKNLKPDKPM QFLGDEETVR KAMEAVAAQG KAKK.