

# Bovine TOP1 Recombinant Protein



RPPB2438

---

## Product Information    Protein Information

**Product SKU:**

RPPB2438

**Protein description:**

Bovine DNA Topoisomerase-I purified from bovine tissues by proprietary chromatographic techniques.

**Accession:**

Q2HPH0

**Synonyms:**

DNA topoisomerase 1, EC 5.99.1.2, DNA topoisomerase I, TOP1, Scl-70.

**Host:**

Bovine tissues.

**Formulation:**

TOP1 is supplied at a concentration of 1mg/1ml in 16mM HEPES buffer pH-7.5, 400mM NaCl and 20% glycerol.

**Purity:**

Greater than 70% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.