

Mouse CAMK2N1 Recombinant Protein



RPPB2448

Product Information Protein Information

Product SKU:

RPPB2448

Accession:

Q6QWF9

Host:

Escherichia Coli.

Protein description:

CAMK2N1 Mouse Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 101 amino acids (1-78 a.a) and having a molecular mass of 10.9kDa. CAMK2N1 is fused to a 23 amino acid His-tag at N-terminus & purified by proprietary chromatographic techniques.

Appearance:

Sterile filtered colorless solution.

Synonyms:

mCaMKII α , calcium/calmodulin-dependent protein kinase II inhibitor α , Camk2n1.

Formulation:

CAMK2N1 protein solution (0.5mg/ml) containing 20mM Tris-HCl (pH8.5), 1mM DTT, 0.1M NaCl and 10% glycerol.

Purity:

Greater than 90.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

Amino Acid Sequence:

MGSSHHHHHH SGLVPRGSH MGSMSEVLPY GDEKLSPLYGD GGDVGGQIFSC RLQDTNFFFG AGQSKRPPKL
GQIGRSKRVV IEDDRIDDLV KTMTDKAPPG V.