

Human CKMBITII Recombinant Protein



RPPB2484

Product Information Protein Information

Product SKU:

RPPB2484

Host:

Pichia Pastoris.

Protein description:

CKMBITII Human Recombinant produced in Pichia Pastoris reacts with polyclonal antibodies to MB Isoenzyme in ELISA.

Appearance:

Sterile Filtered colourless liquid formulation.

Synonyms:

Creatine Kinase MB Isoenzyme Type-II, CKMBITII, CKMBI, CKMB.

Formulation:

Each mg of protein contains 20mM Tris-HCl, pH-6.8, 1mM EDTA and 1mM DTT, 50% (v/v) glycerol.

Stability:

CKMBITII although stable at 15°C for 7 days, should be stored below -18°C. Please prevent freeze-thaw cycles.

Biological Activity:

The enzyme activity measured by kinetic assay at 340nm was 884.6 IU/mg at 37 degrees.