

# Human PPARG (1-477) Recombinant Protein



RPPB2625

---

## Product Information    Protein Information

**Product SKU:**

RPPB2625

**Accession:**

P37231

**Host:**

Escherichia Coli.

**Protein description:**

PPARG Human Recombinant encoding amino acids 1-477 expressed in E.coli, is fused to a GST tag and it is antibody reactive. The PPARG is purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Peroxisome proliferator-activated receptor gamma, PPAR-gamma, PPARG, NR1C3, PPARG1, PPARG2.

**Formulation:**

PPARG at 0.1mg/ml in 50mM Tris-HCl and 10mM L-glutathione (reduced).

**Stability:**

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Avoid multiple freeze-thaw cycles.