

Human AFP Recombinant Protein



RPPB2828

Product Information Protein Information

Product SKU:

RPPB2828

Protein description:

Human alpha-fetoprotein purified from pooled human cord serum.

Accession:

P02771

Appearance:

Sterile Filtered clear solution.

Host:

Human cord serum.

Synonyms:

Alpha-fetoprotein, Alpha-fetoglobulin, Alpha-1-fetoprotein, AFP, FETA, HPAFP.

Formulation:

AFP protein filtered (0.2µm) solution in Tris buffered saline pH 7.5 and less than 0.1% NaN₃.

Stability:

Human Alpha-Fetoprotein should be stored at 2-8°C.