

Bartonella 17kDa Recombinant Protein



RPPB2834

Product Information Protein Information

Product SKU:

RPPB2834

Host:

Escherichia Coli.

Protein description:

Recombinant Bartonella Henselae 17kDa produced in E.Coli is a single, non-glycosylated polypeptide chain having a molecular mass of 18kDa. Bartonella 17kDa is expressed with a -10x His tag at N-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered solution.

Formulation:

Bartonella 17kDa is supplied at a 20mM HEPES buffer pH-7.6, 150mM NaCl and 40% Sucrose (w/v).

Purity:

Greater than 80% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.