

# Bartonella 26kDa Recombinant Protein



RPPB2835

---

## Product Information    Protein Information

**Product SKU:**

RPPB2835

**Host:**

Escherichia Coli.

**Protein description:**

Recombinant Bartonella Henselae 26kDa produced in E.Coli is a single, non-glycosylated polypeptide chain having a molecular mass of 25kDa. Bartonella 26kDa is expressed with a -10x His tag at N-terminus and purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered solution.

**Formulation:**

Bartonella 26kDa is supplied at a 20mM HEPES buffer pH-7.6, 150mM NaCl and 40% Sucrose (w/v).

**Purity:**

Greater than 80% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.