

# Human AMBP Recombinant Protein



RPPB2850

---

## Product Information    Protein Information

**Product SKU:**

RPPB2850

**Accession:**

P02760

**Host:**

Purified from the urine of patients with chronic renal tubular proteinuria.

**Protein description:**

Alpha 1-microglobulin (A1M) is an immunomodulatory protein with a broad spectrum of possible clinical applications and seems a promising marker for evaluation of tubular function.

**Appearance:**

Sterile Filtered White lyophilized (freeze-dried) powder.

**Synonyms:**

Alpha-1 Microglobulin, A1M.

**Formulation:**

Lyophilized from 0.02M NH<sub>4</sub>HCO<sub>3</sub>. May contain traces of buffer salts.

**Purity:**

Greater than 96.0%.

**Solubility:**

Use phosphate buffer, pH>7.0 containing 0.15M NaCl, is recommended.

**Stability:**

Human A1M although stable at room temperature for 3 weeks, should be stored between 2-8°C.