

Mouse Bax Recombinant Protein



RPPB2855

Product Information Protein Information

Product SKU:

RPPB2855

Accession:

Q07813

Host:

Escherichia Coli.

Protein description:

Bax Mouse Recombinant amino acid 1-171 produced in E.Coli is a single, non-glycosylated polypeptide chain. The Mouse Bax is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered colorless solution.

Synonyms:

Apoptosis regulator BAX membrane isoform alpha, Bax, Bcl2-associated X protein.

Formulation:

The 1mg/ml protein contains 10mM Tris-HCl pH-8, 1mM EDTA and 150mM NaCl, 10% glycerol, 0.1% tween & 10mM glutathione.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.