

# Bovine ATF Recombinant Protein



RPPB2893

---

## Product Information    Protein Information

**Product SKU:**

RPPB2893

**Protein description:**

Bovine Apo Transferrin is a glycoprotein of approximately 77kDa.

**Accession:**

Q29443

**Appearance:**

Sterile Filtered off-white lyophilized (freeze-dried) powder.

**Synonyms:**

Serotransferrin, Transferrin, Siderophilin, Beta-1-metal-binding globulin, TF, PRO1557, PRO2086, DKFZp781D0156, Apo Transferrin, ATF.

**Formulation:**

The protein (1mg/ml) was lyophilized with no additives.

**Purity:**

Greater than 95.0% as determined by SDS-PAGE.

**Solubility:**

It is recommended to reconstitute the lyophilized Apo Transferrin in sterile 18MΩ-cm H<sub>2</sub>O not less than 1gr/30ml in 20min. at 20-25C.

**Stability:**

Lyophilized Bovine Apo Transferrin although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Bovine Apo Transferrin should be stored at 4°C between 2-7 days and for future use below -18°C. Please prevent freeze-thaw cycles.