

# Human BP-1 Recombinant Protein



RPPB2904

---

## Product Information    Protein Information

**Product SKU:**

RPPB2904

**Accession:**

Q92988

**Host:**

Escherichia Coli.

**Protein description:**

BP-1 Human Recombinant (aa 103-168) expressed in E.coli, shows a 35 kDa band on SDS-PAGE. The BP-1 is purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Homeobox protein DLX-4, DLX-7, DLX-8, Beta protein 1, BP1, DLX7, DLX8, DLX9, DLX4.

**Formulation:**

BP-1 protein in 50mM Tris-Acetate, pH7.5, 1mM EDTA and 20% Glycerol.

**Stability:**

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.