

Bovine APOH Recombinant Protein



RPPB2964

Product Information Protein Information

Product SKU:

RPPB2964

Accession:

P17690

Host:

Fetal calf serum (FCS).

Protein description:

Bovine Apolipoprotein-H polypeptide is purified from fetal calf serum (FCS) by proprietary protein-chemical techniques, having an Mw of approximately 70kDa.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Beta-2-glycoprotein 1, Apolipoprotein H, Apo-H, Beta-2-glycoprotein I, B2GPI, Beta(2)GPI, APOH.

Formulation:

APOH protein solution is supplied in 16mM Tris-HCl pH-8.0, 280mM NaCl, 0.16mM AEBSF and 20% glycerol.

Purity:

Greater than 70% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.