

# Human B2M Recombinant Protein



RPPB2978

---

## Product Information   Protein Information

**Product SKU:**

RPPB2978

**Accession:**

P61769

**Host:**

Human urine from patients with tubular proteinuria.

**Protein description:**

Human Beta-2 Microglobulin produced in Human urine from patients with tubular proteinuria having a molecular mass of 12KDa.

**Appearance:**

Sterile Filtered White lyophilized (freeze-dried) powder.

**Synonyms:**

Beta-2-microglobulin, B2M.

**Formulation:**

Lyophilized from 0.02M NH<sub>4</sub>HCO<sub>3</sub>.

**Purity:**

Greater than 98.0%.

**Solubility:**

It is recommended to reconstitute the lyophilized HuB2M in phosphate buffer pH > 7.0 containing 0.15M NaCl.

**Stability:**

Human B2M although stable at room temperature for 3 weeks, should be stored between 2-8°C.