

Human CD2 Recombinant Protein



RPPB3015

Product Information Protein Information

Product SKU:

RPPB3015

Accession:

P06729

Host:

Escherichia Coli.

Protein description:

CD2 Human Recombinant (aa 25-209) expressed in E.coli, shows a 46 kDa band on SDS-PAGE. The CD2 is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

T-cell surface antigen CD2, T-cell surface antigen T11/Leu-5, LFA-2, LFA-3 receptor, Erythrocyte receptor, Rosette receptor, CD2 antigen, CD2, T11, SRBC.

Formulation:

CD2 at 100µg/ml in 50mM Tris-HCl pH7.5 and 10Mm L-glutathione (reduced).

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.