

# Human BPI Recombinant Protein



RPPB3045

---

## Product Information    Protein Information

**Product SKU:**

RPPB3045

**Protein description:**

The native protein BPI is isolated from human buffy coats and has a molecular mass of 53kDa.

**Accession:**

P17213

**Synonyms:**

Bactericidal/permeability-increasing protein, CAP 57.

**Formulation:**

BPI (1.15mg/ml) is supplied in 25mM Sodium Acetate pH-4.6.

**Purity:**

Greater than 70% as determined by SDS-PAGE and capillary electrophoresis.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.