

# Human CD34 Recombinant Protein



RPPB3067

---

## Product Information    Protein Information

**Product SKU:**

RPPB3067

**Accession:**

P28906

**Host:**

Escherichia Coli.

**Protein description:**

CD34 Human Recombinant full length expressed in E.coli, CD34 is antibody reactive. The CD34 is purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Hematopoietic progenitor cell antigen CD34, CD34.

**Formulation:**

CD34 in 50mM Tris-HCl, pH7.5, 10mM L-glutathione (reduced).

**Stability:**

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.