

Human CD68 Recombinant Protein



RPPB3086

Product Information Protein Information

Product SKU:

RPPB3086

Accession:

P34810

Host:

Escherichia Coli.

Protein description:

CD68 Human Recombinant full length expressed in E.coli, shows a 38 kDa band on SDS-PAGE. CD68 is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Macrosialin, Gp110, CD68 antigen, CD68, SCARD1, DKFZp686M18236.

Formulation:

CD68 at 100µg/1ml in 50mM Tris-Acetate, pH7.5, 1mM EDTA, 20% Glycerol.

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.