

Human CD72 Recombinant Protein



RPPB3089

Product Information Protein Information

Product SKU:

RPPB3089

Accession:

P21854

Host:

Escherichia Coli.

Protein description:

CD72 Human Recombinant (aa 1-98) expressed in E.coli, shows a 38 kDa band on SDS-PAGE. The CD72 is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

B-cell differentiation antigen CD72, Lyb-2, CD72, LYB2, CD72b.

Formulation:

CD72 at 100µg/ml in 50mM Tris-HCl, pH7.5, 10mM L-glutathione (reduced).

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.