

Human CRP Recombinant Protein



RPPB3411

Product Information Protein Information

Product SKU:

RPPB3411

Accession:

P02741

Host:

Human Plasma.

Protein description:

Human CRP produced in Human plasma having a molecular mass of 114 kDa. It can be used as a marker for inflammation and also used for monitoring and prediction of future events in coronary artery disease.

Appearance:

Sterile Filtered clear solution.

Synonyms:

C-reactive protein, CRP, PTX1, MGC88244, MGC149895.

Formulation:

The protein solution is in 20mM TRIS buffer pH 8.0 containing 0.28M NaCl, 0.09% NaN₃ and 5mM CaCl₂.

Purity:

Greater than 96.0%.

Stability:

Human CRP should be stored at 2-8°C. Do not freeze!