

Human CST3 Recombinant Protein



RPPB3432

Product Information Protein Information

Product SKU:

RPPB3432

Appearance:

Sterile Filtered White lyophilized (freeze-dried) powder.

Accession:

P01034

Synonyms:

Cystatin-C, Cystatin-3, Neuroendocrine basic polypeptide, Gamma-trace, Post-gamma-globulin, CST3, MGC117328.

Host:

Human urine.

Formulation:

Lyophilized from 0.02M NH₄HCO₃. May contain traces of buffer salts.

Purity:

Greater than 96% as determined by SDS PAGE.

Solubility:

Use of a 0.02M NH₄HCO₃ buffer is recommended.

Stability:

Human Cystatin-C although stable at 4°C for 1 week, should be stored at -15°C. Please avoid freeze-thaw cycles.