

Human FABP3 Recombinant Protein



RPPB3495

Product Information Protein Information

Product SKU:

RPPB3495

Accession:

P05413

Host:

Escherichia Coli.

Protein description:

FABP3 Human Recombinant is a non-glycosylated polypeptide chain produced in E.Coli. FABP3 is purified by standard chromatography techniques.

Appearance:

Sterile filtered liquid formulation.

Synonyms:

Fatty acid-binding protein heart, H-FABP, Heart-type fatty acid-binding protein, Muscle fatty acid-binding protein, M-FABP, Mammary-derived growth inhibitor, MDGI, FABP3, FABP11, O-FABP.

Formulation:

50mM phosphate borate buffer pH-8.

Purity:

Greater than 90.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.