

Human FABP3 Recombinant Protein



RPPB3496

Product Information Protein Information

Product SKU:

RPPB3496

Accession:

P05413

Host:

Escherichia Coli.

Protein description:

FABP3 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 133 amino acids and having a molecular mass of 19.1 kDa. FABP3 is fused to a 37 amino acid His tag at N-terminus and purified by standard chromatography techniques.

Appearance:

Sterile Filtered colorless liquid formulation.

Synonyms:

Fatty acid-binding protein heart, H-FABP, Heart-type fatty acid-binding protein, Muscle fatty acid-binding protein, M-FABP, Mammary-derived growth inhibitor, MDGI, FABP3, FABP11, O-FABP.

Formulation:

FABP3 His-Tag is supplied in 20mM Tris-HCl pH 8.0 and 10% glycerol.

Purity:

Greater than 95% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.