

Human FABP4 Recombinant Protein



RPPB3498

Product Information Protein Information

Product SKU:

RPPB3498

Accession:

P15090

Host:

Escherichia Coli.

Protein description:

FABP4 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 132 amino acids and having a molecular mass of 18 kDa. FABP4 is fused to His tag at N-terminus and purified by standard chromatography techniques.

Appearance:

Sterile Filtered colorless liquid formulation.

Synonyms:

Fatty acid-binding protein adipocyte, AFABP, Fatty acid-binding protein 4, Adipocyte lipid-binding protein, ALBP, A-FABP, FABP4.

Formulation:

FABP4 His-Tag is supplied in 20mM Tris HCL pH=8, 0.5mM DTT and 50% glycerol.

Purity:

Greater than 95% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.