

# Human FABP5 Recombinant Protein



RPPB3500

---

## Product Information    Protein Information

**Product SKU:**

RPPB3500

**Accession:**

Q01469

**Host:**

Escherichia Coli.

**Protein description:**

FABP5 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 135 amino acids and having a molecular mass of 19.66kDa. FABP5 is fused to His tag at N-terminus and purified by standard chromatography techniques.

**Appearance:**

Sterile Filtered colorless liquid formulation.

**Synonyms:**

Fatty acid-binding protein epidermal, E-FABP, Fatty acid-binding protein 5, Psoriasis-associated fatty acid-binding protein homolog, PA-FABP, FABP5, EFABP, PAFABP.

**Formulation:**

FABP5 His-Tag is supplied in 20mM Tris HCl pH-8 and 50% glycerol.

**Purity:**

Greater than 95% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.