

# Giardia lamblia Recombinant Protein



RPPB3590

---

## Product Information    Protein Information

**Product SKU:**

RPPB3590

**Host:**

Escherichia Coli.

**Protein description:**

The E.Coli derived recombinant Giardia intestinalis protein contains 266 amino acids (21-280). Giardia lamblia protein is fused to 6xHis tag at C-terminal and purified by proprietary chromatographic techniques.

**Formulation:**

Giardia Intestinalis antigen is formulated at a concentration of 1mg/ml in 1x PBS pH-7.2.

**Purity:**

Giardia lamblia protein is >90% pure as determined by SDS-PAGE.

**Stability:**

Giardia Intestinalis protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.