

Human ICAM-1 Recombinant Protein



RPPB3724

Product Information Protein Information

Product SKU:

RPPB3724

Accession:

P05362

Host:

Escherichia Coli.

Protein description:

ICAM-1 Human Recombinant is expressed in E. coli containing 455 amino acids 28-483 fused to an amino terminal hexahistidine tag. The ICAM-1 is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Intercellular adhesion molecule 1, ICAM-1, Major group rhinovirus receptor, CD54 antigen, ICAM1, BB2, CD54, P3.58.

Formulation:

ICAM-1 is supplied (0.2mg/ml) in PBS, 50% Glycerol, 2mM Imidazole.

Purity:

Greater than 95% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 1-2 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles.