

GCV Recombinant Protein



RPPB3748

Product Information Protein Information

Product SKU:

RPPB3748

Protein description:

Ganciclovir is a white to off-white crystalline powder with a molecular formula of $C_9H_{13}N_5O_4$ and a molecular weight of 255.23.

Appearance:

Sterile Filtered White lyophilized (freeze-dried) powder.

Purity:

Greater than 99.0%.

Solubility:

It is recommended to reconstitute the lyophilized Ganciclovir in sterile 18M Ω -cm H₂O not less than 100 μ g/ml, which can then be further diluted to other aqueous solutions.

Stability:

Lyophilized Ganciclovir although stable at room temperature for 3 weeks, should be stored at 4°C.