

HSA Protein Recombinant Protein



RPPB3890

Product Information Protein Information

Product SKU:

RPPB3890

Protein description:

HSA contains 584 amino acid residues derived from the prototypical HSA sequence.

Accession:

P02768

Appearance:

Sterile Filtered clear yellowish solution.

Host:

Human Serum.

Synonyms:

HSA, ALB, PRO0883, PRO0903, PRO1341, DKFZp779N1935, GIG20, GIG42, PRO1708, PRO2044, PRO2619, PRO2675, UNQ696, SA, HSA.

Formulation:

0.2gr/ml solution containing no additives.

Purity:

Greater than 97.0% as determined by SDS-PAGE.

Stability:

HSA although stable at room temperature for 2 weeks should be stored at 4°C.