

# M.Pneumoniae P1-C Recombinant Protein



RPPB3901

---

## Product Information    Protein Information

**Product SKU:**

RPPB3901

**Host:**

E.Coli

**Protein description:**

The Recombinant Mycoplasma Pneumoniae C-terminal region (P1C) of the P1 protein was expressed in E. coli containing 362 amino acids. The protein is fused to a 6 His Tag. This peptide can be used as an antigen in the diagnosis of M. pneumoniae infection.

**Appearance:**

Sterile filtered colorless solution.

**Formulation:**

Mycoplasma Pneumoniae P1-C is formulated in 1x PBS pH 7.4.

**Purity:**

Greater than 95% as determined by 12% PAGE (Coomassie staining).

**Stability:**

Upon arrival, Store at -20°C. Please prevent freeze-thaw cycles.