

M.Pneumoniae P116 Recombinant Protein



RPPB3902

Product Information Protein Information

Product SKU:

RPPB3902

Host:

E.Coli

Protein description:

The Recombinant Mycoplasma Pneumoniae P116 was expressed in E. coli having an Mw of 116 kDa. P116 is fused to a His-Tag.

Appearance:

Streile filtered colorless solution.

Formulation:

Mycoplasma Pneumoniae P116 is formulated in 1xPBS, pH 7.4 and 25mM arginine.

Purity:

Greater than 95% as determined by 10% PAGE (Coomassie staining).

Stability:

Upon arrival, Store at -20°C. Please prevent freeze-thaw cycles.