

MoeX Recombinant Protein



RPPB3988

Product Information Protein Information

Product SKU:

RPPB3988

Host:

Escherichia Coli.

Protein description:

Recombinant Mycobacterium Tuberculosis MoeX produced in E.Coli is a single, non-glycosylated, polypeptide chain containing having a molecular mass of 38 kDa. The MoeX is fused to a His tag at C-terminus and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Formulation:

PBS.

Purity:

Protein is >95% pure as determined by 10% PAGE (coomassie staining).

Stability:

MoeX Recombinant although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.