

Human LDL Recombinant Protein



RPPB4023

Product Information Protein Information

Product SKU:

RPPB4023

Protein description:

Human Low Density Lipoprotein (LDL) produced in Human plasma.

Host:

Human plasma.

Appearance:

Yellow to orange liquid.

Synonyms:

Low Density Lipoprotein, LDL.

Stability:

Human LDL although stable at 4°C for 1 week, should be stored below -15°C (short term i.e. < 3 months) and below -70°C for long term. Human LDL can be further diluted with saline + 15% sucrose.