

Human MYL1 Recombinant Protein



RPPB4030

Product Information Protein Information

Product SKU:

RPPB4030

Accession:

P05976

Host:

Escherichia Coli.

Protein description:

Recombinant Human Ventricular Myosin Light Chain-1 (MYL1) protein has a molecular mass of 25 kDa and is fused to 7 amino acids at N-terminus. The MYL1 protein was affinity purified using anti MYL1 monoclonal antibody 39-15 column.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Myosin light chain 1 skeletal muscle isoform, MLC1F, A1 catalytic, Alkali myosin light chain 1, MYL1, MLC3F.

Formulation:

Human MYL1 in 10mM Tris -HCl, 1mM EDTA PH 7.5.

Purity:

Greater than 90.0% as determined by SDS-PAGE.

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months. Please prevent freeze-thaw cycles.