

Human Myoglobin Recombinant Protein



RPPB4043

Product Information Protein Information

Product SKU:

RPPB4043

Accession:

P02144

Host:

Escherichia Coli.

Protein description:

Myoglobin Human Recombinant produced in E.Coli is a non-glycosylated polypeptide chain having a molecular mass of 17.67 kDa. The Myoglobin is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered brownish solution.

Synonyms:

Myoglobin, MB, PVALB, MGC13548.

Formulation:

The sterile solution contains phosphate-buffered saline (pH 7.4) and 0.05% NaN₃.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Stability:

Myoglobin should be stored at 4°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please do not freeze.