

Pylori CagA Recombinant Protein



RPPB4160

Product Information Protein Information

Product SKU:

RPPB4160

Host:

Escherichia Coli.

Protein description:

CagA Pylori recombinant antigen is produced in E. coli expressing the H. pylori Cytotoxin-Associated Gene A having the Mw of 34 kDa. ADK is fused to a His-tag at C-terminus & purified by proprietary chromatographic techniques.

Appearance:

Sterile filtered colorless solution.

Formulation:

CagA Pylori recombinant protein is formulated in 1xPBS pH 7.4.

Purity:

Greater than 95% pure as determined by 12% PAGE (Coomassie staining).

Stability:

CagA Pylori although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.