

# S.Enteritidis OmpA Recombinant Protein



RPPB4162

---

## Product Information    Protein Information

**Product SKU:**

RPPB4162

**Host:**

Escherichia Coli.

**Protein description:**

The Recombinant Salmonella Enteritidis Outer Membrane Protein A, E.Coli derived, 330 amino acids, contains the ompA immunodominant regions. The protein is fused to a His tag at C-terminal and purified by standard chromatography techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Outer Membrane Protein-A, OmpA.

**Formulation:**

PBS.

**Purity:**

Protein is >95% pure as determined by 10% PAGE (coomassie staining).

**Stability:**

OmpA S.Enteritidis Recombinant although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.