

Human p53 Recombinant Protein



RPPB4179

Product Information Protein Information

Product SKU:

RPPB4179

Accession:

P04637

Host:

Escherichia Coli.

Protein description:

p53 Human Recombinant full length produced in E.Coli is a non-glycosylated, polypeptide chain having a total Mw of 81kDa. p53 Human Recombinant is fused to GST tag and purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Cellular tumor antigen p53, Tumor suppressor p53, Phosphoprotein p53, Antigen NY-CO-13, TP53, P53, LFS1, TRP53, FLJ92943.

Formulation:

Purified human p53 in 50mM Tris-HCl, pH-7.5 and 10mM L-glutathione (reduced).

Stability:

For long term storage store at -20°C. Avoid freeze/thaw cycles.