

RPPB4499

Product Information Protein Information

Product SKU:

RPPB4499

Protein description:

Prothrombin complex is a double virus-inactivated concentrate, containing equal concentrations of Factors II, VIII, IX, X and Heparin.

Appearance:

Sterile Filtered White lyophilized (freeze-dried) powder.

Synonyms:

PCC, Prothrombin Complex.

Formulation:

Lyophilized Powder from 0.1M sodium chloride and 10mM sodium citrate buffer.

Solubility:

It is recommended to reconstitute the lyophilized Prothrombin Complex in sterile 18M Ω -cm H₂O.

Stability:

Lyophilized Prothrombin-Complex although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Prothrombin-Complex should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.